

Option™ ELITE

Retrievable Vena Cava Filter

KIT COMPONENTS

Introducer Sheath

70cm length 6.5F OD (5F ID) low profile system

Dilator

32mm spaced R/O markers for accurate Vena Cava measurement, dilator side holes arranged in helix formation compatible with power injectors up to 800 psi

Curved Pusher

Allows for controlled and precise over-the-wire filter placement

Filter Cartridge

Universal colour coded cartridge for femoral/jugular placement

Now Available

- High performance delivery sheath with improved kink resistance and pushability
- Radiopaque band for accurate filter placement
- Hemostatic valve

PRODUCT INFORMATION

| Catalogue No | Description | Working Caval Diameter | Introducer Sheath | Length |
|--------------|---------------------------------------|------------------------|-------------------|--------|
| 352506070E | Option™ ELITE Vena Cava Filter System | Up to 32mm | 6.5F OD (5F ID) | 70cm |

EMEA Customer Service

Argon Medical Devices Sarl
Espace Performance
Alphasis Bâtiment J
35760 St Gregoire
France
Tel: +33 (0) 2 99 68 94 95
Fax: +33 (0) 2 99 23 87 21
info-europe@argonmedical.com

Direct Sales Nordics

Argon Medical Devices Danmark A/S
Knud Bro Alle 4
DK-3660 Stenløse
Denmark
Tel: +45 47 10 10 48
Fax: +45 70 40 45 10
customer@argonmedical.com

Argon Medical Devices Danmark A/S
Arenavägen 41 plan 12
SE-121 77 Stockholm
Sweden
Tel: +46 8 50 56 88 70
Fax: +46 8 50 56 88 69
customer@argonmedical.com

Direct Sales France

Argon Medical Devices Sarl
Espace Performance
Alphasis Bâtiment J
35760 St Gregoire
France
Tel: +33 (0) 2 99 68 99 24
Fax: +33 (0) 2 99 23 87 21
service-france@argonmedical.com



Designed By:




© 2013 Argon Medical Devices Inc. All Rights Reserved.
Option™ ELITE is a trademark of Rex Medical.
No reproduction of any portion of this work shall be made without prior permission of the author.
96-6000-30/0713A

Option[™] ELITE

Retrievable Vena Cava Filter

Enhanced for ease of
delivery and retrievability



Proven Safety, Stability
and Structural Integrity

Proven safety and stability required for permanent placement coupled with proven long-term retrievability...

It's your option - Option™ ELITE


Modified Retrieval Apex for ease of snaring

- Single piece construction from nitinol hypotube reduces risk of fracture caused by bending stress (2)
- Transitionless hook allows for ease of snaring during retrieval
- Sturdy design reduces risk of hook deformation during filter retrieval

Enhanced Retention Anchor Pattern for improved delivery and retrievability

- Designed to prevent caudal and cephalad migration
- Retention anchors provide the balance between stability and retrievability
- Curved non-flexible anchors reduce caval wall penetration

The same trusted design enhanced
for ease of delivery and retrievability



Resilient Nitinol Construction

- Resists fracture caused by cyclic stress (1)
- Conforms to caval wall reducing vessel compliance
- Electropolished for increased fracture resistance (1)

Stable Filter Base

- Conforms to anatomy reducing caval penetration (3)
- Provides stability in cavas up to 32mm
- Eliminates caval diameter restrictions
- Aids in centering filter